

## SPEC SHEET

# Drying Tray and Trolley

### INTRODUCTION

The drying tray is used for the final drying of softgel capsules discharged from the softgel capsule tumbler dryer. The Drying Tray ensures capsules remain in the correct shape with the appropriate capsule humidity. Drying trays can be added or removed without affecting the overall running of the cooling device.

**KEYWORDS** Drying tray | trolley | softgel



### TECHNICAL PARAMETERS

Item	Tray (piece)	Trolley
Texture of material	PET	SS304
Dimensions	750×480×53mm	744×476×170mm
Weight	1.65kg	11.6kg

While the information has been carefully compiled to the best of our knowledge, nothing is intended as representation or warranty on our part and no statement herein shall be construed as recommendation to infringe existing patents. In the effort to improve our products, we reserve the right to make changes judged to be necessary.

### CONNECT WITH US

+86 577 6688 6278 | [info@capsulatech.com](mailto:info@capsulatech.com) | [www.captteksolutions.com](http://www.captteksolutions.com)