

SPEC SHEET

Inspection Table

JYT-A MODEL

INTRODUCTION

This machine is the role of the drying complete softgel capsule for effective inspection and selection, to prevent the bubble shaped pill, pill and loading unqualified capsules mixed with the finished product.

KEYWORDS Inspection table | inspection and selection | unqualified softgel capsule | softgel | accessories



ATTRIBUTES

- Waste pill bucket for singling out unqualified capsules
- Lighting lamp for clear inspection

TECHNICAL PARAMETERS

Model	JYT-A
Power Supply	220V 60w
Dimensions	1200×550×1200mm
Weight	65kg

While the information has been carefully compiled to the best of our knowledge, nothing is intended as representation or warranty on our part and no statement herein shall be construed as recommendation to infringe existing patents. In the effort to improve our products, we reserve the right to make changes judged to be necessary.

CONNECT WITH US

+86 577 6688 6278 | info@capsulatech.com | www.captteksolutions.com