

SPEC SHEET

Colloid Mill

JM-130 MODEL

INTRODUCTION

This machine is mill equipment, the main principle is that mill the liquid medicine by gradually reduce the clearance between two grinding wheels under high speed running. The medicine powder can reach super fine 160 mesh after mill.

KEYWORDS Colloid mill | smooth encapsulating | medicine | softgel | accessories



ATTRIBUTES

- Two grinding wheels for fine and even milling
- Cooling device to avoid the temperature rising
- Motor sealed in the machine body to reduce the working noise
- Easy to clean

TECHNICAL PARAMETERS

Model	JM-130
Processing Material Fineness	2~50 μ m
Output (Change for Material Properties)	0.5~2 t/h
Motor Power	380V 50Hz 3P 11kw
Motor Speed	2930 r/min
Rotor Diameter	130mm
Dimensions	Φ 550×1168mm
Weight	320kg

While the information has been carefully compiled to the best of our knowledge, nothing is intended as representation or warranty on our part and no statement herein shall be construed as recommendation to infringe existing patents. In the effort to improve our products, we reserve the right to make changes judged to be necessary.

CONNECT WITH US

+86 577 6688 6278 | info@capsulatech.com | www.captteksolutions.com